1 Einführung

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1 Introduction

1.1 Copyright and Intellectual property

The copyright for these assembly and use instructions is held solely by the manufacturer.
Additionally, all rights are reserved, particularly with regard to patent grant or utility model registration.
Violations of the above-mentioned provisions will result, where applicable, in the obligation to pay compensation for damages.
Intellectual Property – Patent numbers
U.S. Patent No.: 8529484
CA Patent No.: 2692534
European Patent: EP2533734

1.2 Foreword
We would like to thank you for your confidence in TurboMed and the XTERN. This orthosis was designed by TurboMed Orthotics Inc. located in the province of Quebec, Canada. This new-generation orthotic solution was developed specifically for people with foot drop.
The XTERN was developed to offer you unrivaled performance and comfort.
We urge you to read this guide carefully and to address any questions and/or comments to your authorized provider before using or assembling your new orthosis. In order to ensure that the equipment works properly and safely, you must strictly follow all of the instructions provided in this document and observe the specified usage restrictions.
At any time, you have any concerns about the way your device is working or if you have any questions, do not hesitate to contact your authorized provider.
To contact the customer service department of TurboMed Orthotics, please visit: www.TurboMedOrthotics.com
You'll find a complete list of distributors by country on our website: www.TurboMedOrthotics.com

1.3 Use
Intended use: This orthosis is intended SOLELY as an orthotic fitting for the lower limb and has been designed only for drop foot.
It is important to use it in accordance with the indications.

1.4 Safety symbols used in this guide
The symbols below are used in this guide to identify the warnings and useful advice that you should observe. It is very important to read and thoroughly understand them before using our product for the first time.

Warning: this symbol means that there is a potential risk. Failure to heed a warning may result in bodily harm or injury.

Notice: recommendations or instructions should be taken seriously.

Product-related caution: concerns a potential material risk. Failure to heed a product caution may result in significant material damage.

2 Presentation of the product

Item	Description
1	FRONT POSTERIOR SECTION ("Y" SHAPE)
2	FRANT HARNESS ("U" SHAPE)
3	BACK HARNESS ("U" SHAPE)
4	CLASP CLIP
5	CALF SECTION
6	JUNCTION BLOCK

3 Indications and contraindications

3.1 Indications
The orthosis is intended for people who have:
Weakness in the muscles that lift the foot, with little to no spasticity (post: back, hip or knee surgery; traumatic injury to the lower limbs; a stroke; traumatic brain injury; multiple sclerosis; Charcot-Marie-Tooth disease; personal paralysis; neuromuscular atrophy, etc.).
Minor impairment of the muscles that lift the foot.
Very slight deformations of the foot capable of being corrected with the use of an insole and lateral support.
Indications for this product must always be determined by your doctor or prescribing physician.

3.2 Contraindications
The orthosis was not designed for cases of:
-moderate to severe spasticity of the lower limb;
-leg ulcers in the posterior part of the calf;
-severe balance problems;
-moderate to severe foot deformities;
-certain pathologies that require a visit to a physician to verify whether the orthosis is appropriate for the client or not;
-test by the competent authority.

4 The Product
(This product should be used only under the supervision of a qualified physician. MD/DO/DFPM)

4.1 Usage restrictions
Use of the orthosis is not recommended to prevent contractures (shortening of muscles), etc.
Avoid going near a source of high heat, such as a fire, a heating element, embers, flames, etc.
NEVER use an orthosis that belonged to someone else or an orthosis that you believe has been used.

1.2 Warning and general safety instructions (client)

Information:
-please read this document in its entirety before using the product;
-heed the safety instructions in order to avoid injury or damage to the product;
-keep this document in a safe place.

1.3 Warning and general safety instructions (professional)
Information:
-please read this document in its entirety before using the product;
-heed the safety instructions in order to avoid injury or damage to the product;
-show users how to properly use their product and inform them accurately of the safety instructions;
-keep this document in a safe place.

4.4 Caution
(This product is for single patient use only. Use caution when handling the product. Check the product regularly to see if it still works properly. Stop using the product if there are changes in, or loss of, functionality and have it inspected by your authorized health care provider. Observe the legal and insurance regulations related to operating motor vehicles of any type and have your driving license renewed as required. Read the instructions carefully. In general, the physician will determine the daily duration of use and the total period of use. Only qualified health care providers are authorized to perform the initial fitting and first application of the product. Show the patient how to properly handle and care for the product. Inform the patient that it is important to see a physician immediately if unusual changes to the body are observed (e.g., the development of a wound, pain, etc.). Ensure that the orthosis is applied properly and that it fits correctly. Tell the patient to check, prior to each usage, that the product works properly and that it shows no signs of wear or damage. Direct the patient's attention to the fact that they CANNOT continue to use the product if it or one of its components is worn or damaged (ex., cracks, missing part, warping, inappropriate shape of the orthosis, etc.).

5 Application and removal of the product

5.1 Application
1 Install the lace clip to an appropriate shoe.
2 Open the calf band.
3 Install the product onto the orthosis.
4 Put on the shoe and orthosis, pushing back on the orthosis.
5 Close and tie the shoe.
6 Tighten the calf band so as to comfortably secure the top part of the brace to the leg.
Be careful not to over tighten.

5.2 Removal
1 Open the calf band.
2 Untie your shoe.
3 Remove your foot from the shoe.

5.3 Final inspection
Ask the patient to take a few steps to try out the product and to test it on slopes, ramps, while running, and on stairs. Make adjustments, if necessary.

6 Care and cleaning
Care and cleaning should be done on a regular basis. Clean the rigid parts of the product only with mild soap and water or a damp cloth. The so-called "soft" parts, such as the calf band, etc., should be cleaned by hand using mild soap, then air-dried. Avoid any source of direct heat (ex., solar radiation, a radiator, fire, etc.).

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4 The Product
(This product should be used only under the supervision of a qualified physician. MD/DO/DFPM)
4.1 Usage restrictions
Use of the orthosis is not recommended to prevent contractures (shortening of muscles), etc.
Avoid going near a source of high heat, such as a fire, a heating element, embers, flames, etc.
NEVER use an orthosis that belonged to someone else or an orthosis that you believe has been used.

1.2 Warning and general safety instructions (client)
Information:
-please read this document in its entirety before using the product;
-heed the safety instructions in order to avoid injury or damage to the product;
-keep this document in a safe place.

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